


Safety Module

Safety Topic: Piranha Cleaning Solution

Author: Suzanne Szewczyk





- **Equipment** Sulfuric acid, 30% aqueous hydrogen peroxide, glass Erlenmeyer flask, glass waste container, glass graduated cylinder, thick butyl gloves
- **Process** Piranha must be made and used exclusively in a fume hood. Remove all other compounds from the immediate area. Work behind the hood sash/panel. Slowly add 1 portion of hydrogen peroxide to the 3 portions of sulfuric acid with stirring to keep the concentration of peroxide low. The solution may become hot. Carefully add the solution to pre-washed and.dried.glassware to be cleaned. Let sit overnight in the back of a hood away from other materials.
- Specific considerations
 - Do not use plastic containers as they will react with the solution.
 - Piranha generates $O_2(g)$ which can lead to a fire or explosion.
 - Prepare the solution immediately prior to using. Never store piranha and do not put piranha in a <u>closed container or it will explode.</u>
 - Mixing hot piranha with organic compounds such as <u>acetone</u>, <u>isopropanol</u>, <u>nylon</u> <u>or photoresist and may cause an explosion</u>.
 - Neutralize small spills with sodium bicarbonate. A larger spill or a spill outside a fume hood requires EHS assistance for cleanup.
 - Skin or eye exposure to piranha can cause severe burns. The vapor is highly corrosive to mucosal membranes and lungs.
- Waste handling Allow the solution to cool to room temperature and stop generating gasses (~12 hours) before transferring into a clean and dry waste container for disposal. Label the waste as hydrogen peroxide and sulfuric acid.
- For more information see
 - NOCHROMIX® is a safer and more stable cleaning solution alternative (inorganic persulfate solid to be mixed with sulfuric acid) - detailed information on handling and SDS: http://godax.com/msds-directions-for-use
 - o Piranha SOP https://ehs.yale.edu/sites/default/files/files/piranha-sop.pdf
 - O Acid Piranha guidelines https://www.safety.admin.cam.ac.uk/files/hsd176c.pdf